



HSA Approval: MDPR260521F0008

Dengue NS1 & IgG/IgM Rapid Test




Qualitative detection of Dengue NS1 antigen and IgG/IgM antibodies in human capillary whole blood (obtained via finger prick) to aid in the diagnosis of Dengue virus infection.

– WHEN TO TEST

- Sudden high fever ($\geq 38.5^{\circ}\text{C}$)
- Severe headache, pain behind eyes
- Joint & muscle pain
- Recent travel to a dengue-affected area
- Skin rash
- Nausea or vomiting
- Easy bruising or bleeding

– WHEN EACH MARKER APPEARS



– TEST PERFORMANCE

	NS1	IgG	IgM
Sensitivity	97.0%	96.81%	97.87%
Specificity	98.0%	100%	100%